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What is GoutChew Bambi mini ® and when should it be used ?

GoutChew Bambi mini ® is a bioactive oral teat with an ortho-postural aim. Its principle is based on two actions, namely the avoidance of bad habits and the proprioceptive stimulation of the jawbones.

It appears that an increasing number of growing children have difficulty giving up thumb sucking (or other fingers), the lozenge, the bottle or various tissues or favorite stuffed animals. However, when these bad habits persist, they have harmful repercussions on the shape of the jaws, the position of the baby teeth and also the permanent teeth. GoutChew Bambi mini ® has been designed to replace the sucking of a finger (or a foreign body) in order to progressively eliminate the bad habit and allow the resumption of normal growth of the bones of the mouth.

It also appears that the main diet of the youngest children (carefully mixed pre-cooked meals) hardly requires any masticatory effort. The absence of muscular effort and solicitation of the teeth tends to underdevelop the manducatory apparatus. The bones of the mouth are too small, and the space available for the teeth to position themselves correctly on the dental arches is limited. GoutChew Bambi mini ® acts directly on the muscular function of the child's masticatory apparatus by correctly stimulating the growth of the mouth while correcting and guiding the milk teeth and permanent teeth towards a correct positioning on the dental arches.

GoutChew Bambi mini ® should be used when the child :

- Is unable to do without thumb sucking or a foreign object;
- Has a tendency to breathe through the mouth;
- Has a lack of space for the front teeth (overlapping, front teeth);
- Grinds teeth;
- Snoring at night or sleep apnea;
- Has pain related to tooth eruption;
- Is at risk of developing a body posture disorder or prevention.

When should GoutChew Bambi mini ® not be used ?

GoutChew Bambi mini ® is not suitable for children suffering from epilepsy or convulsions.

What precautions should be taken when using GoutChew Bambi mini ® ?

4. Maintain optimal oral hygiene

Brushing your teeth should never be omitted, especially before bedtime, as wearing this device requires reinforcement of oral hygiene.

5. Respect the adaptation time

Wearing a bioactive mouthpiece requires a period of adaptation that can vary from one person to another. This is why it is recommended to wear it gradually, starting with a few minutes a day, until you reach half an hour. If active use causes discomfort in the jaws, it is advisable to use it passively (without exercise) until the situation returns to normal.

6. Adjust the gutter

If the child complains of nausea or the teat is too large, it is advisable to gradually shorten the side branches until there is no more discomfort.

How to use GoutChew Bambi mini® ?

It is advisable to clean the device before use. Wear the bioactive mouthpiece at night and during the day gradually for 5-30 minutes. Perform daily exercises by chewing the mouthpiece with your lips still together. Do not bite the mouthpiece (squeezing does not provide any benefit) but chew it. Avoid talking with the mouthpiece in your mouth.

What side effects can GoutChew Bambi mini ® cause ?

It is possible that the incisors are sensitive at the beginning of the treatment but this disappears quickly. Hypersalivation is also observed and is completely normal. Salivation will gradually decrease and the child will be able to wear the device continuously. There is no need to worry if the child cannot keep it in the mouth all

night, especially if he/she breathes often by mouth. The main thing is the repetition and regularity of the daily exercises to obtain the stimulating effects on growth.

In case of irritation of the palate and hypersensitivity to the teeth, the wearing time should be shortened. In case of an injury to the palate, the orthosis should not be worn until complete healing.

In case of nausea and vomiting, it is possible that the back legs of the orthosis are too long. You then have the option of shortening them with scissors. If the nausea and vomiting still persist, it is likely that the child has a hypersensitive palate. Contact your dentist or speech therapist to see if desensitization is possible.

It may happen that the gums bleed in the initial phase of wearing the splint. In this case, gradual wearing is recommended to ensure that the tissue is strengthened and the bleeding stops spontaneously. The bioactive soother should not cause injury or pain under any circumstances. If it does, ask your dentist to adjust the braces to reduce friction on the mucosa. It is very important to respond to the child and give him rest when discomfort occurs.

What else should we keep in mind ?

1. Maintenance

The bioactive teat should be rinsed with warm water after each use and brushed with toothpaste once a day.

When not in use, it should be stored in its box with the lid open for the first few hours of use to dry.

2. Wear

This therapeutic tool is made of silicone, which loses its elasticity and effectiveness over time, the lifespan of the nipple is approximately 6 months (depending on factors such as activity level, chewing strength, maintenance).

Since the benefits gradually diminish, it is advisable to replace the device after 9 months at the latest.

3. Dental treatment

If the child is undergoing implant, prosthetic or orthodontic treatment, or if the child suffers from periodontitis (loosening of the teeth), ask your dentist for advice before using the splint.

4. Advices

Gradually offer your child harder foods, such as raw vegetables to nibble on or certain fruits to stimulate their chewing system and optimal development.

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