## GoutChew Bambi ®

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## GoutChew Bambi ®

### What is GoutChew Bambi ® and when should it be used ?

GoutChew Bambi ® is a bioactive oral teat with an ortho-postural aim. Its principle is based on two actions, namely the prevention of bad habits and the proprioceptive stimulation of the jaw bones.

An increasing number of growing children seem to have difficulties in giving up thumb sucking (or other fingers), pacifier sucking, bottle sucking, or sucking on various tissues or stuffed animals. If these bad habits persist, they will have a negative effect on the shape of the jaw, the position of the baby teeth and also the permanent teeth. GoutChew Bambi ® has been designed to replace the sucking of a finger (or a foreign body) in order to progressively eliminate this bad habit and to allow a normal growth of the bones of the mouth again. It also appears that the main diet of younger children (carefully mixed ready meals) requires little chewing effort. The lack of muscular effort and the lack of tooth engagement tends to underdevelop the jaw apparatus. Because the bones of the mouth are too small, the space available for the teeth to position themselves correctly in the dental arches is limited. GoutChew Bambi ® acts directly on the muscular function of the child's masticatory apparatus by stimulating the correct growth of the mouth while correcting and guiding the milk teeth and permanent teeth into the correct position in the dental arches.

GoutChew Bambi ® should be used when the child :

- Is unable to do without thumb sucking or a foreign object;
- Has a tendency to breathe through the mouth;
- Has a lack of space for the front teeth (crowding, protruding teeth);
- Teeth grinding;
- Nighttime snoring or sleep apnea;
- Pain from tooth eruption;
- Is at risk of developing a body posture disorder or prevention.

### When should GoutChew Bambi ® not be used ?

GoutChew Bambi ® should not be used in children with epilepsy or seizures.

### What precautions should be taken when using GoutChew Bambi ® ?

#### 4. Maintain optimal oral hygiene

Under no circumstances should you skip brushing your teeth, especially before going to bed, because wearing this device requires reinforced oral hygiene.

#### 5. Respect the adaptation time

Wearing the bioactive mouthpiece requires an adaptation period that can vary from one person to another. It is therefore recommended to wear the mouthpiece progressively, starting with a few minutes a day up to half an hour. If active wearing causes discomfort in the jaw, it is advisable to wear the mouthpiece passively (without moving) until the situation has returned to normal.

#### 6. Adjusting the mouthpiece

If the child complains of nausea or if the nipple is too large, it is recommended to gradually shorten the side branches until the child no longer feels discomfort.

### How to use GoutChew Bambi ® ?

It is recommended to clean the device before using it. Use the bioactive mouthpiece gradually for 5 to 30 minutes in the evening and during the day. Practice daily by chewing the mouthpiece with your lips still together. Do not bite the mouthpiece (squeezing it does not help), but chew it. Avoid talking with the mouthpiece in your mouth.

### What side effects can GoutChew Bambi® cause ?

Incisors may be sensitive at the beginning of treatment, but this disappears quickly. Hypersalivation can also occur, which is normal. The salivation will gradually decrease and the child will be able to wear the appliance continuously. There is no need to worry if the child cannot keep it in the mouth all night, especially if he/she breathes frequently by mouth. The main thing is to repeat the daily exercises regularly in order to obtain the effects favorable to growth.

It is advisable to reduce the duration of wear in case of irritation of the palate and hyper-sensitivity of the teeth. If a lesion of the palate occurs, it is preferable to stop wearing the brace until it has completely healed. In case of nausea and vomiting, it is possible that the posterior branches of the brace are too long. In this case, you can shorten them with scissors. If nausea and vomiting persist, it is likely that the child has a hyper-sensitive palate. Please consult your dentist or speech therapist to see if desensitization is possible. Gums may bleed when the mouthpiece is first worn. A gradual wearing of the mouthpiece is then recommended to ensure a reinforcement of the tissues and the spontaneous stopping of the bleeding. The bioactive teat should not cause any injury or pain. If it does, ask your dental practitioner to adjust the appliance to reduce friction against the mucosa. It is essential to listen to the child and let him/her rest if any discomfort occurs.

#### **What else should we look out for ?**

##### **1. Maintenance**

The bioactive pacifier should be rinsed with warm water after each use and brushed with toothpaste once a day. When not in use, it should be stored in its box with the lid open within the first few hours of use so that it dries.

##### **2. Wear and tear**

This therapeutic tool is made of silicone, which loses its elasticity and effectiveness over time. The lifespan of the orthosis is approximately 6 months (depending on factors such as activity level, masticatory strength, maintenance). As the benefits gradually diminish, it is recommended to replace the appliance after 9 months at the latest.

##### **3. Dental treatments**

If the child is undergoing implant, prosthetic or orthodontic treatment or if he or she suffers from periodontitis (loosening of the teeth), ask your dentist for advice before using the brace.

##### **4. Advice**

Gradually offer your child harder food, such as raw vegetables to nibble on or certain fruits to stimulate their chewing system and optimal development.

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