GoutChew Pingou ®

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GoutChew Pingou ®

What is GoutChew Penguin ® and when should it be used ?

GoutChew Pingou ® is a bioactive oral orthosis with an ortho-postural aim. Its principle is based on two actions: proprioceptive stimulation of the perioral soft tissues and the bones of the mouth, as well as the promotion of balanced muscular work of the jaws, which cascades down to the body posture.

It appears that food in general has become much softer over the decades (cooking, fast food, mixed dishes), requiring almost no masticatory effort. The absence of muscular effort and solicitation of the teeth tends to underdevelop the manducatory apparatus. The bones of the mouth are too small, and the space available for the teeth to position themselves correctly on the dental arches is limited. GoutChew Penguin ® acts directly on the muscular function of the child's masticatory apparatus by correctly stimulating the growth of the mouth while correcting and guiding the permanent teeth to a proper position on the dental arches.

It also appears that the mandibular joint (temporo-mandibular) is a sensor allowing us to situate ourselves in space (see postural system) of as much importance as our eyes or our feet. This joint allows our mandible to move in all directions in space, which makes it the most complex joint in the human body. An imbalance on one side of this double joint can cause the mandible to be misaligned and the body to compensate, resulting in a disturbance of the overall body posture and causing tension. GoutChew Pingou ® relaxes and rebalances the forces acting on the temporomandibular joints, thus freeing the growing child from any compensatory stress that could permanently fix the body.

GoutChew Penguin ® should be used when the child :

- Grinds teeth;
- Interposes his tongue between his teeth when swallowing (persistent infantile swallowing);
- Tends to breathe through the mouth;
- Lack of space for front teeth (overlapping, forward teeth);
- Snoring at night or sleep apnea;
- Feels anxious or stressed;
- Is at risk of developing a body posture disorder or prevention;
- Has pain related to tooth eruption.

When should GoutChew Pingou ® not be used ?

GoutChew Pingou ® is not suitable for children suffering from epilepsy or convulsions. The use of the oral orthosis is also contraindicated :

- If you have significant tooth mobility: if you feel your permanent teeth moving when you apply manual pressure, you should not use the brace;
- If you have an incisor that is significantly longer than the others: difference equal to or greater than 3mm;
- For 3 months following the end of an orthodontic treatment.

What precautions should be taken when using GoutChew Penguin ® ?

4. Maintain optimal oral hygiene

Brushing your teeth should never be omitted, especially before bedtime, since wearing this device requires you to reinforce your oral hygiene.

5. Respect the adaptation time

Wearing the bioactive orthosis requires a period of adaptation that can vary from one person to another.

Therefore, it is recommended to wear it gradually, starting with a few minutes a day, up to half an hour. If active use causes discomfort in the jaws, it is advisable to use it passively (without exercise) until the situation returns to normal.

6. Adjust the gutter

If the child complains of nausea or the brace is too large, it is advisable to gradually shorten the side arms until there is no more discomfort.

How to use GoutChew Penguin ® ?

It is advisable to clean the device before use. Wear the bioactive mouthpiece gradually for 5-30 minutes a day when walking or doing sports. For intense physical efforts requiring mouth breathing, it is advisable to use it just during the warm-up. Perform daily exercises by chewing the mouthpiece with your lips still together. Do not bite

the mouthpiece (clenching does not provide any beneficial effect) but chew it. In order to benefit from a greater muscular relaxation, it is advised to do small liberating exercises after a few minutes of use :

- Turn and tilt the head to each side;
- Turn and tilt the trunk to each side;
- Lean forward, legs straight, hands as close as possible to the ground;
- Squat down and stand up slowly, rolling out the back;
- Repeat these exercises 2-3 times.

After 7 to 10 days, you can start wearing it at night to maximize the benefits, especially if you grind your teeth (bruxism), if you suffer from sleep apnea or if your snoring is annoying.

What side effects can GoutChew Penguin ® cause ?

It is possible that the incisors are sensitive at the beginning of the treatment but this disappears quickly. Hypersalivation is also observed and is completely normal. Salivation will gradually decrease and the child will be able to wear the device continuously. There is no need to worry if the child cannot keep it in the mouth all night, especially if he/she breathes often by mouth. The most important thing is to repeat the daily exercises regularly in order to obtain the stimulating effects on growth.

In case of irritation of the palate and hyper-sensitivity of the teeth, it is advisable to reduce the wearing time. In the case of a lesion on the palate, it is preferable to stop wearing the brace until it has completely healed.

In the case of nausea and vomiting, it is possible that the posterior branches of the orthosis are too long. In this case you can shorten them with scissors. If the nausea and vomiting persists, it is likely that the child has a hyper-sensitive palate. Please refer to your dentist or speech therapist to see if desensitization is possible.

It can happen that the gums bleed at the beginning of the wearing of the mouthpiece. In this case, gradual wear is recommended to ensure that the tissues are strengthened and that the bleeding stops spontaneously. The bioactive brace should not cause any injury or pain. If it does, ask your dental practitioner to adjust the appliance to reduce friction against the mucosa. It is essential to listen to the child and to let him rest if any discomfort appears.

What else should I look out for ?

1. Maintenance

The bioactive oral orthosis must be rinsed with warm water after each use, brushed with toothpaste once a day. When not in use, it should be stored in its box with the lid open within the first few hours of use to allow it to dry.

2. Wear

This therapeutic tool is made of silicone, which loses its elasticity and effectiveness over time. The lifespan of the orthosis is 3 to 6 months (depending on factors such as activity level, masticatory strength, maintenance). As the benefits gradually diminish, it is recommended to replace the appliance after 6 months at the latest.

3. Dental treatments

If you are undergoing implant, prosthetic or orthodontic treatment or suffer from periodontitis (loosening of the teeth), limit active wearing of the brace to a maximum of 10 minutes during the day (can still be worn at night) and ask your dentist for advice before using the brace.

4. Advice

Gradually offer children harder foods, such as raw vegetables to snack on or certain fruits to stimulate their chewing system and optimal development.

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